

Hyalus Disposable Culture Tubes (Borosilicate Glass Tubes)

Description:

Hyalus premium quality glass tubes are produced from borosilicate glass that provides excellent chemical resistance and durability. The tubes are available in Various sizes and feature uniform bottoms and fire-polished rims. They are neatly packaged in boxes of 250 and are shrink-wrapped in plastic film to keep them particle free.

Features and Benefits:

- Made from superior quality borosilicate glass to reduce pH changes and contaminants that are potentially leached from tubes made from soda lime glass.
- Tubes can be used for tissue culture and blood banks.
- Tubes are chemically inert.
- Tubes feature well-rounded bottoms and smooth, fire-polished rims.
- Dimensional specifications are rigidly controlled for reliable, uniform results.

Products:

Catalogue No	Description	Packaging
H- 10075	10 x 75mm, 3mL	250 per box
H- 12075	12 x 75mm, 5mL	250 per box
H- 13100	13 x 100mm, 7mL	250 per box
H- 16100	16 x 100mm, 10mL	250 per box
H- 16125	16 x 125mm, 12mL	250 per box
H- 16150	16 x 150mm, 15mL	250 per box



Technical Specifications:

Material: Type 1, Class B Borosilicate Glass

Autoclavable: Yes

Spin Rate: Max. ~ 4500 rpm

Melting Point (Log n=4): 1460°C (+/-20°C)

Softening Point (Log n=7.6): 1220°C (+/-20°C)


Freezing Point: -50°C

Expansion = $5 \sim 6 \times 10^{-6} \text{ K}^{-1}$

Product Certifications:

- ISO 9001:2008
- Meets ASTM E890 Specification, Type 1 Requirements (Standard Specification for Disposable Glass Culture Tubes)

Note: Any institution or individual who reuses a disposable glass culture tube must bear full responsibility for its safety and effectiveness.



Glass Syringe For Bottle Top Dispenser

Sailent Feartures:

- Fully Customised
- Frosted Glass Barrel and Glass Piston for jerk free Smooth movement
- Borosilicate Glass
- 100% Inert

Sizes:

2.5ml, 5ml, 10ml, 25ml, 60ml, 100ml



Pipette Syringe



Graduated pipette with Suction Piston

- Pipette Complies with IS, ISO / DIN 835
- Pipette are graduated to the tip
- Graduation in white enamel increase the visibility with dark solution

Sizes	Catalogue No	Capacity (ml)	Subdivision (ml)	Tolerance \pm ml	Packing
1ml	GSPS - 01 - 0.01	1	0.01	0.01	1
1ml	GSPS - 01 - 0.1	1	0.1	0.01	1
2ml	GSPS - 02 - 0.02	2	0.02	0.02	1
2ml	GSPS - 02 - 0.01	2	0.1	0.02	1
5ml	GSPS - 05 - 0.05	5	0.05	0.05	1
5ml	GSPS - 0.5 - 0.1	5	0.1	0.05	1
10ml	GSPS - 10 - 0.1	10	0.1	0.1	1
25ml	GSPS - 25 - 0.2	25	0.2	0.2	1



TRÜTH Gas Syringes

Gas Syringes, 50ml & 100ml with Capillary Tube*

* Also Available with 3 way stopcock



ADDRESS:

Top Syring Compound, Western Express Highway, Behind Samarat Hotel,
Next to Aqua Lodha, Pandurang wadi, Mira Rd (E), Thane 401 104, Mumbai, INDIA

Mob: +91 9769175036

Website: sales@top-syringe.com, www.top-syringe.com